

BacSomatic™

Integrated bacteria and somatic cell count



Product	ID	Time	IBC	SCC
Repeatability	498-1-1	13:34		
Milk SCC	497	13:30		5461
Milk IBC SCC	496	13:07	6.70	87.1
Milk IBC SCC	495	13:05	8.20	91.0
Milk IBC SCC	494	13:01	11.9	85.6
Milk IBC SCC	493	12:40	1875	91.9



ANALYTICS BEYOND MEASURE

BacSomatic™ provides rapid hygiene testing of milk. BacSomatic is the first-ever integrated bacteria and somatic cell tester and offers full automation for minimal reagent handling and consistent test results.

Integrated bacteria and somatic cell tester

Automatic and fast alternative to manual assay, providing instant, simultaneous results for individual bacteria and somatic cell count within 9.5 minutes and 1.5 minutes for somatic cell count alone.

More accurate than alternative methods

A fully automated procedure avoids the risk of human error and inconsistency, while ready-to-use reagents in an enclosed bag system, ensure the exact same dosage for every measurement. A sensor in the reagent bag indicates the number of available tests remaining. The unique reagent concept includes control samples for quality assurance and GLP. Online remote monitoring ensures consistent performance and high uptime.

Simple and flexible to use

Built on a long heritage of FOSS analytical technology and backed by global FOSS support, BacSomatic combines robust performance with a modern interface. Very easy-to-use with intuitive touch-screen operation for smarter money-saving decisions on how to use milk.

Sample types

Raw cow's milk, buffalo milk and goat milk

Parameters

Individual Bacteria Count (IBC)
 Somatic Cell Count (SCC)

Technology

Flow cytometry technology, counting individual bacteria cells and somatic cells.

The IBC results can be converted to CFU results. BacSomatic software includes an easy and rapid tool to develop a robust conversion table between IBC and CFU.

Approvals

EURL/Microval
 FDA/NCIMS
 In compliance with ISO/IDF standards

Specifications

Feature	Specification
Analysis time	9.5 minutes for IBC, 1.5 minutes for SCC
Analysis capacity	SCC and IBC: 15 samples/hour IBC only: 15 samples/hour SCC only: 40 samples/hour
Sample intake	IBC & SCC: 7.2 ml IBC: 6.4 ml SCC: 2.5 ml
Sample temperature	2 - 42 °C (35.6-107.6 °F) (preserved) 2 - 4 degrees °C (non preserved)
Sample quality	Raw milk of normal composition and good quality. Unpreserved or preserved with azidol
Carry-over effect	<1% relative, usually <0.5% (uncompensated) IBC and SCC
Working factor	IBC: 100 and SCC: 100
Measuring range	IBC: 5,000-20 mill. IBC/ml (approx. 1,500 to 10 mill CFU/ml) SCC: 0-10 mill. cells/ml
Performance range	IBC: 10,000 - 10 mill. IBC/ml SCC: 100,000 - 1.5 mill. cells/ml
Accuracy	IBC: Typical $S_{y,x} < 0.25$ log units from SPC (standard plate count) SCC: < 10% relative mean difference form Direct Microscopic SCC
Reference or anchor method	Standard Plate Count (SPC) EN-ISO 4833-1:2013
Dimensions (w x h x d)	400 x 400 x 400 mm
Weight	25 kg
Space requirements (w x d)	Approx. 600 x 630 mm
Power supply	110 - 240 V AC, 50/60 Hz
Power consumption	max 190 VA
Water supply for preparation	Purified water (<5 $\mu\text{S}/\text{cm}^3$) The instrument does not require a direct connection to water tap but for chemical preparations ultra pure de-mineralized water must be available
Waste	Max 1 l
Ambient temperature	Ambient temperature 5 - 35 °C (41 - 95 °F)
Relative humidity	<93% RH
Instrument management	
Networking software	FossManager™

Repeatability		
Range (IBC/ μl)	S_r (log-units)	Typical S_r (log units)
10 - 50	0.07	0.06
51 - 200	0.05	0.04
>200	0.04	0.02
Entire range	0.05	
Reproducibility (between instruments)		
Range (IBC/ μl)	S_R (log-units)	Typical S_R (log units)
10 - 50	0.11	0.08
51 - 200	0.07	0.06
>200	0.06	0.04

Repeatability	
Range (SCC/ μl)	CV, %
100 - 299	6
300 - 499	4
>500	3
Reproducibility	
Range (SCC/ μl)	S_{Rr} %
100 - 299	7
300 - 499	5
>500	4

000 «Диаэм»

Москва
ул. Магаданская, д. 7, к. 3 ■ тел./факс: (495) 745-0508 ■ sales@dia-m.ru

www.dia-m.ru

С.-Петербург
+7 (812) 372-6040
spb@dia-m.ru

Новосибирск
+7 (383) 328-0048
nsk@dia-m.ru

Воронеж
+7 (473) 232-4412
vrn@dia-m.ru

Йошкар-Ола
+7 (927) 880-3676
nba@dia-m.ru

Красноярск
+7 (923) 303-0152
krsk@dia-m.ru

Казань
+7 (843) 210-2080
kazan@dia-m.ru

Ростов-на-Дону
+7 (863) 303-5500
rnd@dia-m.ru

Екатеринбург
+7 (912) 658-7606
ekb@dia-m.ru

Кемерово
+7 (923) 158-6753
kemerovo@dia-m.ru

Армения
+7 (094) 01-0173
armenia@dia-m.ru

