

## BagPipet®

Fixed-volume pipet & straw tips

Ref 251 091 (BagPipet®)



Standards:



Manufactured under:



BagPipet® is an autoclavable pipet with fixed volume (0.1 mL / 0.9 mL / 1 mL). Used with BagTips® a single use rigid straw tip, available in 2 sizes (Regular 19 cm / Jumbo 24 cm), BagPipet® intakes any product, even very viscous products: milk, cream, yoghurt, cosmetics...The long BagTips® reaches the bottom of any bag.

### TECHNICAL SPECIFICATIONS

- Fixed volume pipet
- 3 settings: 0.1 mL / 0.9 mL / 1 mL
- Calibrated according to the GLP "Good Laboratory Practice"
- Autoclavable
- Ejection of disposable straw
- Easy calibration
- Light and ergonomic
- Also available: variable volume button from 10 to 1000 µL (Ref 251 080)
- Box (w x d x h): 30 x 10 x 6 cm, gross weight: 250 g
- 1 year warranty
- In compliance with: ISO 7218, ISO 6887, ISO 4833-2 and FDA BAM (Bacteriological Analytical Manual)
- Designed and made in France

### DELIVERED WITH

- 1 user's manual
- 1 pouch BagTips® Regular
- 1 fixed volume button 1000 µL (on the pipet)
- 1 holder for key
- 2 fixed volume buttons (900 µL, 100 µL)
- 2 O-Rings
- 1 Allen wrench
- Oil

### ACCESSORIES

- Ref. 252 019 **BagTips® Regular** : Straw tip (5 to 12 kGy gamma ray treated) for BagPipet.  
Nominal volume: 1.1 mL - Length: 190 mm – Quantity per box: 1000  
Box (w x d x h): 27,5 x 23 x 11 cm, gross weight: 1,2 kg
- Ref. 252 024 **BagTips® Jumbo** : Straw tip (5 to 12 kGy gamma ray treated) for BagPipet.

# interscience

Nominal volume: 2,2 mL – Length: 240 mm - Quantity per box: 1000  
Box (w x d x h): 27,5 x 23 x 11 cm, gross weight: 1,5 kg

**000 «Диаэм»**

**Москва**

ул. Магаданская, д. 7, к. 3 ■ тел./факс: (495) 745-0508 ■ [sales@dia-m.ru](mailto:sales@dia-m.ru)

[www.dia-m.ru](http://www.dia-m.ru)

**С.-Петербург**  
+7 (812) 372-6040  
[spb@dia-m.ru](mailto:spb@dia-m.ru)

**Новосибирск**  
+7 (383) 328-0048  
[nsk@dia-m.ru](mailto:nsk@dia-m.ru)

**Воронеж**  
+7 (473) 232-4412  
[vrn@dia-m.ru](mailto:vrn@dia-m.ru)

**Йошкар-Ола**  
+7 (927) 880-3676  
[nba@dia-m.ru](mailto:nba@dia-m.ru)

**Красноярск**  
+7 (923) 303-0152  
[krsk@dia-m.ru](mailto:krsk@dia-m.ru)

**Казань**  
+7 (843) 210-2080  
[kazan@dia-m.ru](mailto:kazan@dia-m.ru)

**Ростов-на-Дону**  
+7 (863) 303-5500  
[rnd@dia-m.ru](mailto:rnd@dia-m.ru)

**Екатеринбург**  
+7 (912) 658-7606  
[ekb@dia-m.ru](mailto:ekb@dia-m.ru)

**Кемерово**  
+7 (923) 158-6753  
[kemerovo@dia-m.ru](mailto:kemerovo@dia-m.ru)

**Армения**  
+7 (094) 01-0173  
[armenia@dia-m.ru](mailto:armenia@dia-m.ru)

