

BagRack®

Storage rack for bags

Ref 221 010 (80/100 mL) - **Ref 221 040** (400 mL) - **Ref 221 350** (3500 mL - 4 compartments)
Ref 221 352 (3500 mL - 2 compartments)



Manufactured under:



BagRack® enables the clean and methodical storage of bags.

TECHNICAL SPECIFICATIONS

- All stainless steel
- Autoclavable
- Stackable
- Available for 80/100/400/2000/3500 mL bags
- Designed and made in France

BagRack® 80/100

Ref: 221 010

- For 11 bags 80/100 mL
- Dimensions (w x d x h): 13 x 27.5 x 15 cm
- Box (w x d x h): 28.5 x 14.5 x 16 cm, gross weight: 0.8 kg

BagRack® 400

Ref: 221 040

- For 10 bags 400 mL
- Dimensions (w x d x h): 20.5 x 35.5 x 21 cm
- Box (w x d x h): 37 x 22 x 22 cm, gross weight: 1.4 kg

BagRack® 3500 - 4 compartments

Ref: 221 350

- For 4 bags of 2000/3500 mL
- Dimensions (w x d x h): 43 x 37 x 37 cm
- Box (w x d x h): 50.5 x 41.5 x 42.5 cm, gross weight: 3 kg

BagRack® 3500 - 2 compartments

Ref: 221 352

- For 2 bags of 2000/3500 mL
- Dimensions (w x d x h): 35.5 x 20.5 x 20 cm
- Box (w x d x h): 39 x 29 x 26 cm, gross weight: 1.3 kg

OPTIONS / ACCESSORIES

- Ref. 221 042 **DripTray®**: Stainless steel drip tray for **BagRack® 400** and **BagRack® 3500** – 2 compartments
Box (w x d x h): 36 x 22 x 3 cm, gross weight: 1 kg

000 «Диаэм»

Москва

ул. Магаданская, д. 7, к. 3 ■ тел./факс: (495) 745-0508 ■ sales@dia-m.ru

www.dia-m.ru

С.-Петербург
+7 (812) 372-6040
spb@dia-m.ru

Новосибирск
+7 (383) 328-0048
nsk@dia-m.ru

Воронеж
+7 (473) 232-4412
vrn@dia-m.ru

Йошкар-Ола
+7 (927) 880-3676
nba@dia-m.ru

Красноярск
+7 (923) 303-0152
krsk@dia-m.ru

Казань
+7 (843) 210-2080
kazan@dia-m.ru

Ростов-на-Дону
+7 (863) 303-5500
rnd@dia-m.ru

Екатеринбург
+7 (912) 658-7606
ekb@dia-m.ru

Кемерово
+7 (923) 158-6753
kemerovo@dia-m.ru

Армения
+7 (094) 01-0173
armenia@dia-m.ru

