

2021

www.herenz.de

ΔΙΑ•Μ
современная диагностика

www.dia-m.ru
диагн-он-лайн



Medical and Laboratory Supplies Wash bottles



Safety wash bottles with imprint

Narrow-mouth, PE-LD (PP for Acetone), transparent, spray tube from PP, with coloured VENT-CAP which supports easy visual identification of the media. More safety due to the durable safety imprint in accordance to Directive (EC) No. 1272/2008 (GHS), Material name in German, English, French and Spanish, Chemical formula, CAS-No., hazard pictogram, signal word, risk phrases (H phrases), safety phrases (P phrases), as well as NFPA-Code.

Item	Volume	GL	Height	Ø	VENT-CAP	MOQ
			mm	mm	colour	
• Imprint „Acetone“						
1083270	500 ml	25	180	74	red	12 pcs/carton
1083271	1.000 ml	32	215	92	red	12 pcs/carton
• Imprint „Ethanol“						
1083272	500 ml	25	180	74	orange	12 pcs/carton
1083273	1.000 ml	32	221	92	red	12 pcs/carton
• Imprint „Isopropanol“						
1083274	500 ml	25	180	74	yellow	12 pcs/carton
1083275	1.000 ml	32	221	92	red	12 pcs/carton



Wash bottles - with imprint „Distilled Water“

Narrow-mouth bottles, PE-LD, transparent, screw cap and spray tube from PP, imprint in German, English, French and Spanish.

Item	Volume	Thread	Height	Ø	MOQ
		(GL)		mm	
1083210	250 ml	25	135 mm	65	12 pcs/carton
1083211	500 ml	25	180 mm	74	12 pcs/carton
1083212	1.000 ml	32	221 mm	92	12 pcs/carton



Wash-bottles, PE-LD

Without impression, round, translucent, narrow neck, Tube and screw-cap PP made.

Item	Volume	Thread	Ø	Height	PU = MOQ
1083200	50 ml	GL 18	37 mm	85 mm	24 pcs/carton
1083201	100 ml	GL 18	43 mm	114 mm	24 pcs/carton
1083202	250 ml	GL 25	58 mm	135 mm	12 pcs/carton
1083203	500 ml	GL 25	74 mm	180 mm	12 pcs/carton
1083204	1000 ml	GL 32	92 mm	221 mm	12 pcs/carton



Wash-bottles, coloured

Narrow-mouth bottles, made from PE-LD. Available in four different colours. Screw cap and spray tube made from PP. Thread GL 25.

Item	Volume	Height	Ø	PU = MOQ
	ml	mm	mm	
• red				
1083215	250	135	65	5 pcs/carton
1083216	500	180	74	5 pcs/carton
• green				
1083217	250	135	65	5 pcs/carton
1083218	500	180	74	5 pcs/carton
• yellow				
1083219	250	135	65	5 pcs/carton
1083220	500	180	74	5 pcs/carton
• blue				
1083221	250	135	65	5 pcs/carton
1083222	500	180	74	5 pcs/carton
• Set: 1 item each, 500 ml, blue, yellow, red und green				
1083223	500	180	74	1 Set/carton



000 «Диаэм»

Москва
ул. Магаданская, д. 7, к. 3 ■ тел./факс: (495) 745-0508 ■ sales@dia-m.ru

www.dia-m.ru

С.-Петербург
+7 (812) 372-6040
spb@dia-m.ru

Новосибирск
+7 (383) 328-0048
nsk@dia-m.ru

Воронеж
+7 (473) 232-4412
vrm@dia-m.ru

Йошкар-Ола
+7 (927) 880-3676
nba@dia-m.ru

Красноярск
+7 (923) 303-0152
krsk@dia-m.ru

Казань
+7 (843) 210-2080
kazan@dia-m.ru

Ростов-на-Дону
+7 (863) 303-5500
rnd@dia-m.ru

Екатеринбург
+7 (912) 658-7606
ekb@dia-m.ru

Кемерово
+7 (912) 158-6753
kemerovo@dia-m.ru

Армения
+7 (094) 01-0173
armenia@dia-m.ru

