

# CULTURA Mini-Incubator

ДИАМ  
современная лаборатория

www.dia-m.ru  
заказ on-line

## Multi-functional:

for the incubation of

- 18 dip-slides or
- 12 Petri-dishes or
- blood culture bottles or
- identification systems (any kind)

## Ease of use:

- fixed-temperature 35°C
- variable control (range: 25-40°C)
- visual control guaranteed (clear Plexiglass-door)
- easy to clean

## Reliable:

- high build quality
- space saving
- safe handling (VDE/TÜV);  
CE and GS-certified acc. to international regulations)
- maintenance free

## Complete:

- all in one pack (Mini-Incubator, Multirack, thermometer)

## Key-stone:

- central principle of a Microbiology-laboratory

## Technical data

External dimensions:	W x H x D = 310 x 155 x 168
Internal dimensions:	W x H x D = 220 x 120 x 150
Weight (net):	1.1 kg
Consumption:	26 W/220 V 26 W/110 V
Temperature range:	25-45 °C (fix: 35 °C)
Temperature accuracy:	± 1 °C

## Ordering Information

Product	Merck Cat. No.	Pack size
CULTURA Mini-Incubator (26 W/220 V)	1.13311.0001	1 Incubator 1 Multirack 1 Thermometer
CULTURA Mini-Incubator (26 W/110 V)	1.15533.0001	1 Incubator 1 Multirack 1 Thermometer
MULTIRACK (spare-rack) for CULTURA Mini-Incubator	1.13312.0001	



CULTURA Mini-Incubator

000 «Диаэм»

Москва  
ул. Магаданская, д. 7, к. 3 ■ тел./факс: (495) 745-0508 ■ sales@dia-m.ru

www.dia-m.ru

С.-Петербург  
+7 (812) 372-6040  
spb@dia-m.ru

Новосибирск  
+7 (383) 328-0048  
nsk@dia-m.ru

Воронеж  
+7 (473) 232-4412  
vrn@dia-m.ru

Йошкар-Ола  
+7 (927) 880-3676  
nba@dia-m.ru

Красноярск  
+7 (923) 303-0152  
krsk@dia-m.ru

Казань  
+7 (843) 210-2080  
kazan@dia-m.ru

Ростов-на-Дону  
+7 (863) 303-5500  
rnd@dia-m.ru

Екатеринбург  
+7 (912) 658-7606  
ekb@dia-m.ru

Кемерово  
+7 (923) 158-6753  
kemerovo@dia-m.ru

Армения  
+7 (094) 01-0173  
armenia@dia-m.ru

