




schuett **petriturn**



Turntables for Petri dishes & Inoculating hooks

schuett **biotec**.de

Easy & Uniform

Comfortable



schuett petritum-E with foot-switch and turntable-attachment for Petri dishes up to Ø 150 mm (optional)



Inoculating hook made of stainless steel (Drigalski-hook)



Inoculating hook made of glass (Drigalski-hook)

Petri dishes turntables

For uniform plating of liquids on agar-surfaces in Petri dishes. Easy to handle, tireless.

schuett petritum-E
Electrically driven Petri dish turntable

With electrically drive and sensor. Body made of stainless-steel. Handling of the turntable contactless by sensor. For Petri dishes up to Ø 100 mm (up to Ø 150 mm optional).

Continuously adjustable, constant speed of 10-100 rpm and continuously adjustable operating time of 3 s to 2 min or permanent running.

Optional connection of a foot-switch.

Turntables for Petri dishes of other diameters available on request.

schuett petritum-M
Manual operation

Powerful and long-running due to high mass of turntable.

Useable on both sides, that means two work spaces in one device by easiest turn-over. For Petri dishes with Ø 100 and Ø 150 mm diameter.



schuett petritum-M with inoculating hook (optional)

Technical data

petritum-E

Dimensions (Ø x h)	160 x 70 mm
Weight	approx. 1 kg
Wall power supply	12 V DC 5 VA
Speed continuously, adjustable	up to 10-100 rpm
Operating time, continuously adjustable	3 s - 2 min / permanent running

petritum-M

Dimensions (Ø x h)	160 x 45 mm
Weight	approx. 2.3 kg

Ordering Information

schuett petritum-E, electrically driven, 230 V	Cat.-No. 3.361 252
schuett petritum-E, electrically driven, 115 V	3.361 262
schuett petritum-M, manual operation	3.361 102

Accessories

Foot-switch for schuett petritum-E	3.361 302
Turntable-attachment for Petri dishes up to Ø 150 mm	3.361 402

Inoculating hook

Inoculating hook made of stainless steel (Drigalski-hook), 180 mm length	3.362 152
Inoculating hook made of glass (Drigalski-hook), 145 mm length	3.362 162

000 «Диаэм»

Москва
ул. Магаданская, д. 7, к. 3 ■ тел./факс: (495) 745-0508 ■ sales@dia-m.ru

www.dia-m.ru

С.-Петербург
+7 (812) 372-6040
spb@dia-m.ru

Новосибирск
+7(383) 328-0048
nsk@dia-m.ru

Воронеж
+7 (473) 232-4412
vrn@dia-m.ru

Йошкар-Ола
+7 (927) 880-3676
nba@dia-m.ru

Красноярск
+7(923) 303-0152
krsk@dia-m.ru

Казань
+7(843) 210-2080
kazan@dia-m.ru

Ростов-на-Дону
+7 (863) 303-5500
rnd@dia-m.ru

Екатеринбург
+7 (912) 658-7606
ekb@dia-m.ru

Кемерово
+7 (923) 158-6753
kemerovo@dia-m.ru

Армения
+7 (094) 01-0173
armenia@dia-m.ru

