

Lab jack with metal knob



Lab jack with autoclavable metal knob

Properties:

- Special feature: Autoclavable through metal knob!!!
- All models are manufactured according to DIN 12897.
- Improved scissor technology with reduced tolerances
- Ergonomic knob, rounded and user-friendly
- Scissor-type design made of 18/10 steel

Regarding the operating load, we distinguish between maximum dynamic operating load (kg max dyn) and maximum static operating load (kg max stat) of the lab jack. These standards are also in line with DIN 12897.

Материал	A (Длина)	B (Ширина)	C мин	C макс	Kg max dyn	Kg max stat	Резьба	Арт. Номер
Нержавею- щая сталь 18/10	160	130	60	275	7	30	M8 x 1,25	11123
Нержавею- щая сталь 18/10	200	200	60	275	7	30	M8 x 1,25	11133
Нержавею- щая сталь 18/10	240	240	60	275	7	30	M8 x 1,25	11143