

ELMI CM-50 Micro Centrifuge



Product Description

CM-50 is one of the most compact and powerful micro-centrifuges available today. Its patented air cooling system and magnetic rotor cap guarantee exceptionally low sample heating, far surpassing the industry's most stringent standards. This mighty mini-centrifuge spins down samples at 15,000 RPM, yet runs on only 24 volts, and fits easily into the smallest spaces. Its high-performance motor with maintenance-free drive ensures outstanding durability and long life in your lab.

Features

- 15,000 RPM
- 15,300 RCF
- Low sample heating – max 10°C while running at max speed
- Rapid spin mode
- Secure lid lock – lid unlocks automatically when rotation is complete
- Digital timer and sound notification
- 5 rotor brake settings
- Rotor imbalance sensor
- Patented magnetic rotor cap
- RPM/RCF optional indicator
- Extremely quiet at top speed

Specifications

Centrifugation Speed	1,000-15,000 RPM
Maximum Centrifugal Force	15,294 RCF
Timer	from 0.1 up to 99 min
Modes/Dimensions (L x W x H)	200 x 180 x 145 mm
Power Supply Adapter	24 VDC, 5A
Weight	3.1 kg



Accessories



0.5ml test tubes adapter for Rotor #50.01 and #50.02

Hermetic rotor adapter for 0.5mL test tubes.
Compatible with ELMI CM-50 and CM-50M.



0.2ml test tubes adapter for Rotor #50.01 and #50.02

Hermetic rotor adapter for 0.2mL test tubes.
Compatible with ELMI CM-50 and CM-50M.



Fugamix Series Rotors 50.01

Universal high speed rotor with patent pending cover cap. Insures optimum air flow and incredibly low sample heating.
Total Number of Tubes: 12
Max Tube Capacity (mL): 0.2 - 2ml
Top Speed (RPM): 15,000



Fugamix Series Rotors 50.03

Universal PCR strip rotor for standard PCR strips.
Total Number of Tubes: 2 PCR Strips
Max Tube Capacity (mL): 0.2ml
Top Speed (RPM): 12,500

000 «Диаэм»

Москва
ул. Магаданская, д. 7, к. 3 ■ тел./факс: (495) 745-0508 ■ sales@dia-m.ru

www.dia-m.ru

С.-Петербург
+7 (812) 372-6040
spb@dia-m.ru

Новосибирск
+7(383) 328-0048
nsk@dia-m.ru

Воронеж
+7 (473) 232-4412
vrn@dia-m.ru

Йошкар-Ола
+7 (927) 880-3676
nba@dia-m.ru

Красноярск
+7(923) 303-0152
krsk@dia-m.ru

Казань
+7(843) 210-2080
kazan@dia-m.ru

Ростов-на-Дону
+7 (863) 303-5500
rnd@dia-m.ru

Екатеринбург
+7 (912) 658-7606
ekb@dia-m.ru

Кемерово
+7 (923) 158-6753
kemerovo@dia-m.ru

Армения
+7 (094) 01-0173
armenia@dia-m.ru

