

Microcentrifuge Tubes

Corning offers two styles of microcentrifuge tubes: traditional snap cap tubes for quick access or screw cap tubes for greater sealing security.

Snap Cap Polypropylene Microcentrifuge Tubes

- ▶ Costar® microcentrifuge tubes are certified RNase-/DNase-free
- ▶ Supplied nonsterile and are autoclavable
- ▶ External graduations and frosted writing spot for easy sample identification
- ▶ Positive seal design allows for repeated opening and closing
- ▶ Flat cap surface for convenient labeling
- ▶ Withstands a maximum RCF of 17,000 x g
- ▶ Costar low binding microcentrifuge tubes feature a bonded polymer technology that reduces protein and nucleic acid binding, resulting in better sample recovery



Microcentrifuge Tubes

Snap Cap Polypropylene Microcentrifuge Tubes Ordering Information

Cat. No.	Volume (mL)	Color	Qty/Pk	Qty/Cs
<i>Snap Cap Microcentrifuge Tubes</i>				
3208	0.65	Natural	500	1,000
3209*	0.65	Rainbow*	200	1,000
3620	1.7	Natural	500	500
3621	1.7	Natural	500	5,000
3622*	1.7	Rainbow*	100	500
3213	2.0	Natural	500	1,000
<i>Low Binding Snap Cap Microcentrifuge Tubes Ordering Information</i>				
3206	0.65	Natural	500	500
3207	1.7	Natural	250	250


*Rainbow pack includes one bag each of blue, green, yellow, red, and orange tubes.

Screw Cap Polypropylene Microcentrifuge Tubes


- ▶ Corning® polypropylene microcentrifuge tubes feature screw caps that provide a tight secure seal
- ▶ Choice of attached cap with silicone O-ring or unattached rim seal cap
- ▶ All tubes have a large white marking spot.
- ▶ Withstands a maximum RCF of 13,000 x g
- ▶ Sterile



Microcentrifuge Tubes



▶ Attached loop cap allows for optimum one-handed convenience. Silicone O-ring gasket provides a snug seal, safeguarding samples against leakage.



▶ Easy-to-use unattached rim seal cap design twists on or off in a single turn.

Screw Cap Polypropylene Microcentrifuge Tubes Ordering Information

Cat. No.	Volume (mL)	Cap Style	O-ring	Self Standing	Qty/Cs
430909	1.5	Attached	Yes	No	500
430915	2.0	Attached	Yes	Yes	500
430917	2.0	Unattached	No	Yes	500

000 «Диаэм»

Москва
ул. Магаданская, д. 7, к. 3 ■ тел./факс: (495) 745-0508 ■ sales@dia-m.ru

www.dia-m.ru

С.-Петербург
+7 (812) 372-6040
spb@dia-m.ru

Новосибирск
+7 (383) 328-0048
nsk@dia-m.ru

Воронеж
+7 (473) 232-4412
vrn@dia-m.ru

Йошкар-Ола
+7 (927) 880-3676
nba@dia-m.ru

Красноярск
+7 (923) 303-0152
krsk@dia-m.ru

Казань
+7 (843) 210-2080
kazan@dia-m.ru

Ростов-на-Дону
+7 (863) 303-5500
rnd@dia-m.ru

Екатеринбург
+7 (912) 658-7606
ekb@dia-m.ru

Кемерово
+7 (923) 158-6753
kemerovo@dia-m.ru

Армения
+7 (094) 01-0173
armenia@dia-m.ru

