



DURAN  
WHEATON  
KIMBLE

Excellence in your hands

# WHEATON®

## AUTOMATIC CRIMPING TOOL BATTERY POWERED



### Battery Powered Crimping Tools

Our battery powered crimpers and decappers fit standard laboratory vials and seals with aluminum or two-part aluminum/ steel caps. The electronically adjustable cycle is accurate and reproducible and can be adjusted with the control pad on top of the unit. The tools are ergonomically designed to reduce strain and arm injury associated with many manual crimping tools.

Rechargeable: The tools can be recharged in one to two hours with the universal power supply included. The power supply includes a plug set for operation in most countries.

| Product Highlights   |  |
|--|--|
| • Powerful motor with a fast cycle   | • Intuitive LCD screen makes adjusting the crimp easy                          |
| • Long lasting battery with Lithium Ion technology - up to 800 crimps per charge | • Advanced charging system including warning icon when recharging is necessary |
| • Can be operated while charging   | • Optional base and mounting kit   |
| • For flip-off caps, available in 13 mm and 20 mm sizes*                         |  |

\*Not for use with Flip-up caps

# AUTOMATIC CRIMPING TOOL — BATTERY POWERED

## Battery Powered Crimping Tool

| Cat. No. | Description                      | Size (mm) | Qty |
|----------|----------------------------------|-----------|-----|
| W225808  | Battery Powered Crimper          | 8         | 1   |
| W225811  | Battery Powered Crimper          | 11        | 1   |
| W225812  | Battery Powered Decapper         | 11        | 1   |
| W225813  | Battery Powered Crimper          | 13        | 1   |
| W225814  | Battery Powered Decapper         | 13        | 1   |
| W225820  | Battery Powered Crimper          | 20        | 1   |
| W225821  | Battery Powered Decapper         | 20        | 1   |
| W225815  | Battery Powered Flip Cap Crimper | 13        | 1   |
| W225822  | Battery Powered Flip Cap Crimper | 20        | 1   |

## Automatic Crimping Tool Stand

| Cat. No. | Description         | Qty/Case |
|----------|---------------------|----------|
| W225701  | Crimping Tool Stand | 1        |

## Stoppers and Seals

| Cat. No.        | Description   | Qty/Case |
|-----------------|---|----------|
| <b>Stoppers</b> |   |          |
| W224100-400     | 13 mm Ultra Pure straight plug stopper                    | 1000     |
| W224100-405     | 20 mm Ultra Pure straight plug stopper                    | 1000     |
| <b>Seals</b>    |   |          |
| 224177-01       | 13 mm open top unlined aluminum seal, natural             | 1000     |
| 224182-01       | 13 mm center disc tear-out unlined aluminum seal, natural | 1000     |
| 224192-01       | 13 mm complete tear-off unlined aluminum seal, natural    | 1000     |
| 224202          | 13 mm flip cap unlined aluminum seal, red                 | 1000     |
| W224207         | 13 mm flip cap tear-off unlined aluminum seal, red        | 1000     |
| 224178-01       | 20 mm open top unlined aluminum seal, natural             | 1000     |
| 224183-01       | 20 mm center disc tear-out unlined aluminum seal, natural | 1000     |
| 224191          | 20 mm solid top unlined aluminum seal, natural            | 1000     |
| 224193-01       | 20 mm complete tear-off unlined aluminum seal, natural    | 1000     |
| 224203          | 20 mm flip cap unlined aluminum seal, red                 | 1000     |
| 224208          | 20 mm flip cap tear-off unlined aluminum seal, red        | 1000     |



**000 «Диаэм»**

Москва  
ул. Магаданская, д. 7, к. 3 ■ тел./факс: (495) 745-0508 ■ sales@dia-m.ru

[www.dia-m.ru](http://www.dia-m.ru)

**С.-Петербург**  
+7 (812) 372-6040  
spb@dia-m.ru

**Новосибирск**  
+7 (383) 328-0048  
nsk@dia-m.ru

**Воронеж**  
+7 (473) 232-4412  
vrn@dia-m.ru

**Йошкар-Ола**  
+7 (927) 880-3676  
nba@dia-m.ru

**Красноярск**  
+7 (923) 303-0152  
krsk@dia-m.ru

**Казань**  
+7 (843) 210-2080  
kazan@dia-m.ru

**Ростов-на-Дону**  
+7 (863) 303-5500  
rnd@dia-m.ru

**Екатеринбург**  
+7 (912) 658-7606  
ekb@dia-m.ru

**Кемерово**  
+7 (923) 158-6753  
kemerovo@dia-m.ru

**Армения**  
+7 (094) 01-0173  
armenia@dia-m.ru

