



Gerber Instruments

NEW – Centrifuge Gerber Compact III



ДИА•М
современная лаборатория

www.dia-m.ru
заказ on-line

GERBER
INSTRUMENTS
Since 1892 

Centrifuge Gerber Compact III

Gerber Instruments, since its founding in 1892, has developed and produced analytical devices and instruments for quality control in the dairy industry. For the determining the fat content in milk and milk products as well as the solubility of milk powders centrifuges are the basic equipment in every laboratory.

With our new multipurpose centrifuge Compact III these methods can be carried out:

- Fat determination according to Gerber and Babcock
- Solubility determination according to Gerber and ADPI

The Gerber centrifuge Compact III with high safety standards and reliability and these product features:

- Strong sheet steel housing with an acid-resistant paint work
- Cover lock to prevent opening of the centrifuge during operation
- Unbalance shut-off through uneven filling or sudden glass break
- Motor with elastic bearing for smooth run
- Increased precision thanks to controlled speed



Centrifuge Compact III with universal rotor and Butyrometer / solubility sleeves



Universal rotor with 8 solubility inserts ADPI

Specifications

Control	Microprocessor control
Timer	0-99 min
Heating	600 W adjustable up to 65 °C
Gerber rotor	with 8 or 12 Butyrometer sleeves
Universal rotor	without inserts (order separately)
Voltage	115 / 230 VAC, 50 / 60 Hz
Dimensions (W × D × H)	430 × 535 × 285 mm
Weight	23 kg without rotor and sleeves
Warranty	2 years

Order No.	Type
03.4100	Centrifuge Compact III / U incl. Universal rotor with 8 places without inserts (order separately)
03.4101	Butyrometer / solubility inserts Gerber 8 sleeves made of brass and 8 reduction rings made of stainless steel
03.4102	Solubility inserts ADPI 8 sleeves and 8 reduction rings made of stainless steel
03.4103	Babcock inserts 8 sleeves made of stainless steel
03.4200	Centrifuge Compact III / 8 incl. Gerber rotor and 8 Butyrometer / solubility sleeves made of brass
03.4300	Centrifuge Compact III / 12 incl. Gerber rotor and 12 Butyrometer / solubility sleeves made of brass

000 «Диаэм»

Москва
ул. Магаданская, д. 7, к. 3 ■ тел./факс: (495) 745-0508 ■ sales@dia-m.ru

www.dia-m.ru

С.-Петербург
+7 (812) 372-6040
spb@dia-m.ru

Новосибирск
+7 (383) 328-0048
nsk@dia-m.ru

Воронеж
+7 (473) 232-4412
vrn@dia-m.ru

Йошкар-Ола
+7 (927) 880-3676
nba@dia-m.ru

Красноярск
+7 (923) 303-0152
krsk@dia-m.ru

Казань
+7 (843) 210-2080
kazan@dia-m.ru

Ростов-на-Дону
+7 (863) 303-5500
rnd@dia-m.ru

Екатеринбург
+7 (912) 658-7606
ekb@dia-m.ru

Кемерово
+7 (923) 158-6753
kemerovo@dia-m.ru

Армения
+7 (094) 01-0173
armenia@dia-m.ru

