

BagSeal®

Thermal sealing unit

Ref 261 000



Manufactured under:



BagSeal® allows a wide and tight sealing of all kind of bags (with/without filter). Perfect closure is guaranteed.

TECHNICAL SPECIFICATIONS

- All stainless steel except handles, welding blades (anodized aluminum), bar joint (chromated galvanized steel)
- Large sealing: 250 x 5 mm
- Adjustable power heating
- LED light code at the end of the sealing (welding process duration indicator: welding + cooling time)

- For bags 80/100/400/2000 mL
- Pressure: 22 Newton/cm²

- Voltage and frequency: 100-127 V~ / 220-240 V~ 50-60 Hz
- Max power: 1600 W
- International protection marking: 40

- Dimensions (w x d x h): 35.9 x 27.6 x 32.2 cm, weight: 10.8 kg
- Box (w x d x h): 49 x 40 x 49 cm, gross weight: 13.6 kg

- Warranty 1 year

- Manufactured under CE, RoHS, WEEE
- Designed and made in France

OPTIONS / ACCESSORIES

- Ref 261 500 Maintenance kit for BagSeal®: 2 sealing threads + 2 fuses + 2 Teflon strips

000 «Диаэм»

Москва

ул. Магаданская, д. 7, к. 3 ■ тел./факс: (495) 745-0508 ■ sales@dia-m.ru

www.dia-m.ru

С.-Петербург
+7 (812) 372-6040
spb@dia-m.ru

Новосибирск
+7 (383) 328-0048
nsk@dia-m.ru

Воронеж
+7 (473) 232-4412
vrn@dia-m.ru

Йошкар-Ола
+7 (927) 880-3676
nba@dia-m.ru

Красноярск
+7 (923) 303-0152
krsk@dia-m.ru

Казань
+7 (843) 210-2080
kazan@dia-m.ru

Ростов-на-Дону
+7 (863) 303-5500
rnd@dia-m.ru

Екатеринбург
+7 (912) 658-7606
ekb@dia-m.ru

Кемерово
+7 (923) 158-6753
kemerovo@dia-m.ru

Армения
+7 (094) 01-0173
armenia@dia-m.ru

