



UV TUBES


Lamps 6W

Model	Wavelength (nm)	Length (mm)
T-40.L	365	1197
T-40.M	312	1197
T-40.C	254	1197
T-30.L	365	892
T-30.C	254	892
T-20.L	365	586
T-20.C	254	586
T-15.L	365	435
T-15.M	312	435
T-15.C	254	435
T-8.L	365	288
T-8.M	312	288
T-8.C	254	288
T-6.L	365	211
T-6.C	254	211

UV TUBES

Spectral Energy Distribution

365nm




UV-A

Peak wavelength : 365nm

6, 8, 15, 20, 30, 40 watt available

G5 & G13 Cap

312nm




UV-B

Peak wavelength : 312nm

8, 15, 40 watt available

G5 & G13 Cap

254nm



UV-C

Peak wavelength : 254nm

6, 8, 15, 20, 30, 40 watt available

G5 & G13 Cap

000 «Диаэм»

Москва

ул. Магаданская, д. 7, к. 3 ■ тел./факс: (495) 745-0508 ■ sales@dia-m.ru

www.dia-m.ru

С.-Петербург
+7 (812) 372-6040
spb@dia-m.ru

Новосибирск
+7(383) 328-0048
nsk@dia-m.ru

Воронеж
+7 (473) 232-4412
vrn@dia-m.ru

Йошкар-Ола
+7 (927) 880-3676
nba@dia-m.ru

Красноярск
+7(923) 303-0152
krsk@dia-m.ru

Казань
+7(843) 210-2080
kazan@dia-m.ru

Ростов-на-Дону
+7 (863) 303-5500
rnd@dia-m.ru

Екатеринбург
+7 (912) 658-7606
ekb@dia-m.ru

Кемерово
+7 (923) 158-6753
kemerovo@dia-m.ru

Армения
+7 (094) 01-0173
armenia@dia-m.ru

