



DURAN
WHEATON
KIMBLE

Excellence in your hands

WHEATON®

AUTOMATIC CRIMPING TOOL

HIGH PERFORMANCE



High Performance Crimping Tools

Our High Performance Crimping Tool is designed with inter-changeable jaws, so that you only need one tool and purchase the size crimper and decappers needed for your laboratory. The jaw sets are available as crimper and decappers which fit standard lab vials seals with aluminum or two-part aluminum/steel caps. This crimping tool has the ability to store separate programs for different caps and cap sizes. The tools are ergonomically designed to reduce strain and arm injury associated with many manual crimping tools.

Product Highlights	
• Our fastest and most powerful crimping tool	• Strong enough for all steel and magnetic caps
• Designed with external power source and cord (no battery)	• Uses interchangeable jaw sets
• Optional crimper stand available.	• Intuitive LCD screen makes adjusting the crimp easy
• For flip-off caps, available in 13 mm and 20 mm sizes*	

*Not for use with Flip-up caps

AUTOMATIC CRIMPING TOOL — HIGH PERFORMANCE

High Performance Crimping Tool

Cat. No.	Jaw(s) Included	Size (mm)	Qty/Case
W225830-[]*	—	—	1
W225831-[]*	Crimper	20	1
W225832-[]*	Crimper & Decapper	20	1
W225833-[]*	Flip Cap Crimper & Decapper	20	1

*Refer to adjacent plug chart to order your proper plug style.

High Performance Jaw Sets

Cat. No.	Jaw Type	Size (mm)	Qty/Case
W225741	Crimper Jaw Set	8	1
W225751	Crimper Jaw Set	11	1
W225752	Decapper Jaw Set	11	1
W225761	Crimper Jaw Set	13	1
W225762	Decapper Jaw Set	13	1
W225763	Flip Cap Crimper Jaw Set	13	1
W225771	Crimper Jaw Set	20	1
W225772	Decapper Jaw Set	20	1
W225773	Flip Cap Crimper Jaw Set	20	1

Automatic Crimping Tool Stand

Cat. No.	Description	Qty/Case
W225701	Crimping Tool Stand	1

Stoppers and Seals

Cat. No.	Description	Qty/Case
Stoppers		
W224100-400	13 mm Ultra Pure straight plug stopper	1000
W224100-405	20 mm Ultra Pure straight plug stopper	1000
Seals		
224177-01	13 mm open top unlined aluminum seal, natural	1000
224182-01	13 mm center disc tear-out unlined aluminum seal, natural	1000
224192-01	13 mm complete tear-off unlined aluminum seal, natural	1000
224202	13 mm flip cap unlined aluminum seal, red	1000
W224207	13 mm flip cap tear-off unlined aluminum seal, red	1000
224178-01	20 mm open top unlined aluminum seal, natural	1000
224183-01	20 mm center disc tear-out unlined aluminum seal, natural	1000
224191	20 mm solid top unlined aluminum seal, natural	1000
224193-01	20 mm complete tear-off unlined aluminum seal, natural	1000
224203	20 mm flip cap unlined aluminum seal, red	1000
224208	20 mm flip cap tear-off unlined aluminum seal, red	1000

	Plug Style
[A]	North American Cord, 120V
[B]	Japan Cord, 100V
[C]	Europe Cord, 230V
[D]	UK Cord, 230V
[F]	Australia / China, 240V
[G]	Italy / Chile, 230V
[I]	India Cord, 230V



000 «Диаэм»

Москва
ул. Магаданская, д. 7, к. 3 ■ тел./факс: (495) 745-0508 ■ sales@dia-m.ru

www.dia-m.ru

С.-Петербург
+7 (812) 372-6040
spb@dia-m.ru

Новосибирск
+7(383) 328-0048
nsk@dia-m.ru

Воронеж
+7 (473) 232-4412
vrn@dia-m.ru

Йошкар-Ола
+7 (927) 880-3676
nba@dia-m.ru

Красноярск
+7(923) 303-0152
krsk@dia-m.ru

Казань
+7(843) 210-2080
kazan@dia-m.ru

Ростов-на-Дону
+7 (863) 303-5500
rnd@dia-m.ru

Екатеринбург
+7 (912) 658-7606
ekb@dia-m.ru

Кемерово
+7 (923) 158-6753
kemerovo@dia-m.ru

Армения
+7 (094) 01-0173
armenia@dia-m.ru

