



STORAGE VIALS

Be Simple, Be Professional

- > A full range storage vials, both OEM&ODM are available.
- > Available as both clear and amber glass vials, screw caps with PTFE/silicone septa, PE lined and PE/Alu Foil.
- > Extensively in environmental, food, and pharmaceutical applications, suitable for the packaging of various pharmaceutical intermediates, high value-added chemicals, chemical reagents, biological reagents, cosmetics, essences and oils, etc.
- > Suitable for long-term storage and transportation of products, has excellent sealing performance.

Volume	Size(mm)	Finish	Color	Qty/Pack	Without Caps		
					Vials ONLY	PP Cap PTFE/Silicone Septa	PP Cap PE Septa
Part No.	Part No.	Part No.					
2ml	12x32	9-425	Clear	100	C0000006	C0000153	C0000159
2ml	12x32	9-425	Amber	100	C0000041	C0000156	C0000159
3ml	16x33	13-425	Clear	100	C0000042	C0000183	C0000184
3ml	16x33	13-425	Amber	100	C0000043	C0000187	C0000188
4ml	15x45	13-425	Clear	100	C0000024	C0000183	C0000184
4ml	15x45	13-425	Amber	100	C0000044	C0000187	C0000188
5ml	18x40	15-425	Clear	100	C0000045	C0000191	C0000195
5ml	18x40	15-425	Amber	100	C0000929	C0000193	C0000194
6ml	16x50	13-425	Clear	144	C0001530	C0000183	C0000184
7ml	22x40	18-400	Clear	100	C0000048	C0000220	C0000222
7ml	22x40	18-400	Amber	100	C0000049	C0000221	C0000223
8ml	17x60	15-425	Clear	100	C0000051	C0000191	C0000195
8ml	17x60	15-425	Amber	100	C0000052	C0000193	C0000194
10ml	22x52	18-400	Clear	100	C0000053	C0000220	C0000222
10ml	22x52	18-400	Amber	100	C0000054	C0000221	C0000223
12ml	19x65	15-425	Clear	100	C0001370	C0000191	C0000195
12ml	19x65	15-425	Amber	100	C0001371	C0000193	C0000194
12ml	19x65	15-425	Amber	144	C0001280	C0000193	C0000194
15ml	21x70	18-400	Clear	144	C0001410	C0000220	C0000222
15ml	21x70	18-400	Amber	144	C0001638	C0000221	C0000223
16ml	22x72.5	18-400	Clear	100	C0000055	C0000220	C0000222
16ml	22x72.5	18-400	Amber	100	C0000056	C0000221	C0000223
20ml	27.5x57	24-400	Clear	100	C0000057	C0000216	C0000218
20ml	27.5x57	24-400	Amber	100	C0000058	C0000215	C0000217
22ml	23x85	20-400	Clear	100	C0001222	C0000227	C0000228
22ml	23x85	20-400	Amber	100	C0001223	C0000989	C0000228
30ml	27.5x84	24-400	Clear	100	C0000059	C0000216	C0000218
30ml	27.5x84	24-400	Amber	100	C0000060	C0000215	C0000217
40ml	27.5x95	24-400	Clear	100	C0000061	C0000216	C0000218
40ml	27.5x95	24-400	Amber	100	C0000062	C0000215	C0000217
60ml	27.5x140	24-400	Clear	100	C0000065	C0000216	C0000218
60ml	27.5x140	24-400	Amber	100	C0000066	C0000215	C0000217

www.alwsci.com
+86-571-88837380
marketing@alwsci.com



COMBINATION PACK

Volume	Size(mm)	Finish	Vial Color	PP Cap Color	Septa Material	Qty/Pack	Part No.
2ml	12x32	9-425	Clear	Black	PTFE/Silicone	100	C0000996
2ml	12x32	9-425	Amber	Black	PTFE/Silicone	100	C0001353
3ml	16x33	13-425	Clear	Black	PE	100	C0001629
4ml	15x45	13-425	Clear	Black	PTFE/Silicone	100	C0000997
4ml	15x45	13-425	Amber	Black	PTFE/Silicone	100	C0001227
4ml	15x45	13-425	Clear	Red	PTFE/Silicone	144	C0001316
4ml	15x45	13-425	Clear	Black	PTFE/Silicone	144	C0001394
4ml	15x45	13-425	Amber	Red	PTFE/Silicone	144	C0001317
5ml	18x40	15-425	Clear	Black	PE	100	C0001551
5ml	18x40	15-425	Amber	Black	PTFE/Silicone	100	C0001791
7ml	22x40	18-400	Clear	Black	PE	100	C0001630
8ml	17x60	15-425	Clear	Black	PTFE/Silicone	144	C0001318
8ml	17x60	15-425	Amber	White	PTFE/Silicone	144	C0001319
10ml	22x52	18-400	Clear	White	PTFE/Silicone	100	C0000967
10ml	22x52	18-400	Amber	White	PTFE/Silicone	100	C0000968
10ml	22x52	18-400	Clear	Black	PTFE/Silicone	100	C0001483
10ml	22x52	18-400	Amber	Black	PE	100	C0001681
12ml	19x65	15-425	Clear	Black	PTFE/Silicone	100	C0001678
20ml	27.5x57	24-400	Amber	White	PTFE/Silicone	100	C0001417
20ml	27.5x57	24-400	Amber	Black	PE	100	C0001680
20ml	27.5x57	24-400	Clear	Black	PTFE/Silicone	144	C0001320
20ml	27.5x57	24-400	Amber	White	PTFE/Silicone	144	C0001321
20ml	27.5x57	24-400	Clear	White	PTFE/Silicone	144	C0001328
20ml	27.5x58	22-400	Clear	White	PTFE/Silicone	100	C0000465
20ml	27.5x58	22-400	Clear	White	PE	100	C0000091
20ml	27.5x58	22-400	Clear	White	Aluminum Foil	100	C0000802
20ml	27.5x58	22-400	Amber	White	PTFE/Silicone	100	C0000404
20ml	27.5x58	22-400	Amber	White	PE	100	C0001097
30ml	27.5x84	24-400	Clear	Black	PTFE/Silicone	100	C0001757
30ml	27.5x84	24-400	Clear	Black	PE	100	C0001631
40ml	27.5x95	24-400	Clear	White	PTFE/Silicone	100	C0001002
40ml	27.5x95	24-400	Clear	Black	PTFE/Silicone	100	C0001143



ЛАБОРАТОРНОЕ ОБОРУДОВАНИЕ ■ ПЛАСТИК ■ СТЕКЛО ■ РЕАКТИВЫ ■ НАБОРЫ

www.dia-m.ru
заказ онлайн