

Reagent Reservoirs



4870 and 4871
50 mL Reagent Reservoir



4872 and 4873
100 mL Reagent Reservoir

Costar® Reagent Reservoirs are ideal for repetitively filling multichannel pipettors

- ▶ Manufactured from modified polystyrene
- ▶ Sterile
- ▶ Disposable

Reagent Reservoirs Ordering Information

Cat No.	Volume (mL)	Color	Qty/Pk	Qty/Cs
4870	50	Natural	5/bag	200
4871	50	Natural	1/bag	100
4872	100	White	5/bag	200
4873	100	White	1/bag	100

Transtar-96® Well Liquid Transfer System



Transtar-96 System

- ▶ The Costar Transtar-96 System is a portable, autoclavable liquid handling device for use with 96 well plates
- ▶ A sterile 96 tip disposable cartridge, which loads into the Transtar system, enables liquids to be aspirated, transferred and dispensed over a volume range of 25 to 200 μ L in 5 μ L increments
- ▶ The Transtar-96 System is ideal for changing cell culture media and screening monoclonal antibodies
- ▶ Transtar-96 System accuracy is rated at $\pm 5\%$ at all volume levels

Transtar-96 Well Liquid Transfer System Ordering Information

Cat. No.	Description	Sterile	Qty/Pk	Qty/Cs
7605	Transtar-96, adjustable-volume pipettor	N/A	1	1
7606	Transtar elevator	N/A	1	1
7610	Transtar disposable cartridges	Yes	1	24
4876	Transtar disposable reservoir liner, open	Yes	1	100
4877	Transtar disposable reservoir liner, 12-channel	Yes	1	100
4878	Transtar disposable reservoir liner, 8-channel	Yes	1	100



Transtar Disposable Cartridge

000 «Диаэм»

Москва

ул. Магаданская, д. 7, к. 3 ■ тел./факс: (495) 745-0508 ■ sales@dia-m.ru

www.dia-m.ru

С.-Петербург
+7 (812) 372-6040
spb@dia-m.ru

Новосибирск
+7 (383) 328-0048
nsk@dia-m.ru

Воронеж
+7 (473) 232-4412
vrn@dia-m.ru

Йошкар-Ола
+7 (927) 880-3676
nba@dia-m.ru

Красноярск
+7 (923) 303-0152
krsk@dia-m.ru

Казань
+7 (843) 210-2080
kazan@dia-m.ru

Ростов-на-Дону
+7 (863) 303-5500
rnd@dia-m.ru

Екатеринбург
+7 (912) 658-7606
ekb@dia-m.ru

Кемерово
+7 (923) 158-6753
kemerovo@dia-m.ru

Армения
+7 (094) 01-0173
armenia@dia-m.ru

